

How to test Otiom in 4 steps

1 Setup Otiom

1. Ensure that you have a new Otiom Starter Kit
2. Download and Login to the Otiom-app. Your login information was sent by email
3. Follow the installation instructions in the Otiom-app

Find our installation guide here:
www.otiom.com/install

Important!
Otiom is ready once you have received a text message informing you that it is activated and ready for use!

2 Settings


1. Make sure that the **Otiom Home base** is located in the chosen test area/office
2. For testing purpose use safety level 6

Find our Personal setting guide here:
www.otiom.com/install

Check that you turned on the **Bluetooth** and **Location** connection and have **internet access** on your smartphone.

3 Test

1. Hang the **Otiom-tag** around your neck and go for a walk (minimum 10 min walk for a good experience)
2. After you leave your **Home base** you will get an alarm
3. Within 5 minutes outside you be able to follow your **Otiom-tag** in the app
4. When you get home, Otiom will automatically cancel the alarm



4 Extra features

1. Locate Indoor
 - Use the "Distance to Otiom-user" function at the top of the map
 - This feature is always active
2. Safetyzone
 - After testing Safetylevel 6 - you may want to try testing Safetylevel 3,4 or 5
 - Remember to add your personal safetyzone in the Otiom-app

